

# AZscript™ Reverse Transcriptase

AZscript Reverse Transcriptase (RT) is an RNA-dependent DNA polymerase used to synthesize complementary DNA (cDNA) from an RNA template.

While originating from Moloney Murine Leukemia Virus (M-MLV), AZscript RT has been engineered for increased thermostability and reduced RNase H activity.

In cDNA synthesis a higher reaction temperature can reduce RNA secondary structures and improve efficiency and specificity of the reaction. Furthermore, a reduction in RNase H activity increases performance of the cDNA synthesis, especially for longer transcripts.

AZscript RT shows excellent performance and sensitivity in RT-qPCR (Fig 1).



Increased  
thermostability



Reduced RNase H  
activity




Fig 1. AZscript RT performance in RT-qPCR

**AZscript shows excellent performance and sensitivity in RT-qPCR.** A serial dilution of MS2 Armored RNA (Asuragen, 52000) was heat lysed (75°C, 5min) and used as template for cDNA synthesis using gene-specific primers, according to manufacturer's instructions. A volume equivalent to 2  $\mu$ l of the cDNA was amplified in a 25  $\mu$ l reaction using TaqPath™ qPCR Master Mix, CG. Data showing amplification plots (A) and standard curve (B) for the dilution series corresponding to 100 000 down to 10 copies of RNA, in triplicate measurements.

## Specifications

|                            |                                            |
|----------------------------|--------------------------------------------|
| Source                     | Recombinantly produced in <i>E. coli</i> . |
| Size                       | 79.5 kDa                                   |
| Temperature optimum (>80%) | 38 - 47°C                                  |
| Inactivation               | 70°C, 15 minutes                           |

## Storage and stability

The enzyme is stable at -20°C for up to 12 months in the supplied storage buffer. Additional data on stability of AZscript RT is available on request.

## Quality control

|                                    |                                                                                                                                                                                                                         |
|------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>dsDNA endonuclease activity</b> | 200 U AZscript RT was incubated with a supercoiled plasmid (1 µg) for 4 hours at 37°C. Agarose gel electrophoresis did not reveal any transformation of closed circular DNA to nicked DNA.                              |
| <b>ssDNA endonuclease activity</b> | 200 U AZscript RT was incubated with M13 ssDNA (0.5 µg) for 4 hours at 37°C. Agarose gel electrophoresis did not reveal any visible signs of ssDNA degradation.                                                         |
| <b>Exonuclease activity</b>        | 200 U AZscript RT was incubated with either 3H-dATP labelled ds or ssDNA (0.5 µg, 500 bp) for 4 hours at 37°C. Acid soluble radioactivity from labelled DNA was not significantly over blank test for either substrate. |
| <b>RNase activity</b>              | 200 U AZscript RT was incubated with a 2 kb RNA transcript (1 µg) for 4 hours at 37°C. Agarose gel electrophoresis did not reveal any visible signs of RNA degradation.                                                 |

## Ordering information

|                     | Article no. | Pack Size* | Concentration |
|---------------------|-------------|------------|---------------|
| <b>AZscript™ RT</b> | 72100-201   | 10 000 U   | 200 U/µl      |
|                     | 72100-100   | Custom     | Custom        |

\*One unit is defined as the amount of enzyme needed to incorporate 1 nmol of dTTP into acid-precipitable material in 10 min at 37°C using poly(A)●oligo(dT)15 as template-primer. The enzyme is assayed in 37.5 mM Tris-HCl pH 8.3 @ 25°C, 40 mM KCl, 6 mM MgCl<sub>2</sub>, 0.75 mg/ml Poly(A), 0.03 mM oligo(dT)15, 10 mM DTT, 0.5 mM dTTP, 0.025 mCi/ml 3H-dTTP, 0.1 mg/ml BSA, 0.02% Triton X-100.

# Your OEM partner to deliver novel solutions for genomics and proteomics.

## Quality

ArcticZymes is dedicated to the quality of its products and is certified according to ISO 13485. ArcticZymes offers the convenience of providing standard bulk enzymes as off the shelf products. In addition, ArcticZymes offers enzymes in customized formats. For additional information, please contact us.

## Contact information

ArcticZymes Technologies ASA  
Sykehusvegen 23, 9019 Tromsø, Norway

**T** +47 7764 8900

**E** [contact@arcticzymes.com](mailto:contact@arcticzymes.com)

**I** [www.arcticzymes.com](http://www.arcticzymes.com)

## Additional Information

We are pleased to provide additional data and information relating to AZscript reverse transcriptase on request. For more information about our enzymes and services, please visit our website [www.arcticzymes.com](http://www.arcticzymes.com).

**Disclaimer:** ArcticZymes Technologies ASA's, including its subsidiaries, ("ArcticZymes") products are intended for the further use of in manufacturing new products or research only. Certain applications of ArcticZymes' products may require licenses from third parties. It is the express duty of any receiver of our products to acquire such licenses, if necessary. To the maximum extent allowed by law, ArcticZymes will not be accountable or liable for any damages, whether direct, indirect, incidental, or consequential, in connection with or arising from this document, including the use of it. ArcticZymes products may be covered by pending or issued patents, designs or design applications and/or trademarks or trade- mark applications or any other registered or unregistered Intellectual Property Right. For the avoidance of doubt, the General Conditions of ArcticZymes shall apply with respect to any and all purchases and use of the ArcticZymes products. Version 1.1 • Feb 2025