



ArcticZymes High-Quality Kits

For contamination
control

www.arcticzymes.com

ArcticZymes High-Quality Kits

For contamination control

ArcticZymes offers novel and high-quality enzymes with unique combination of low temperature activity, as well as heat lability, to facilitate new possibilities in molecular biology applications.



**PCR Decontamination
Kit**



**Heat&Run[®] gDNA
Removal Kit**

PCR Decontamination Kit

- ✓ Reduces background and improves target detection
- ✓ No negative effect on PCR sensitivity
- ✓ Efficient for ordinary PCR and probe-based qPCR mixes



"For effective removal of DNA contamination from PCR master mixes"

PCR is a sensitive method for detecting presence of DNA in both research and diagnostics. Taq polymerases are frequently contaminated with *E. coli* DNA, which might cause reduced sensitivity and false positives when small amounts of bacterial DNA are targeted. Other sources of contamination might be dNTPs, buffer components and primers/probes, as well as DNA introduced during handling.

Description

ArcticZymes' PCR Decontamination Kit offers an easy and affordable solution for removal of contaminating DNA. It utilizes a double-stranded specific dsDNase to remove contaminating DNA. The specificity of dsDNase allows decontamination to occur with primers and probe present. In addition, it is efficient for some SYBR based qPCR mixes.


The dsDNase is irreversibly inactivated by heating to 60°C in presence of DTT, thus ensuring that any template added after inactivation remains safe from digestion.

Specifications

Source	dsDNase is recombinantly expressed in yeast
Contents	dsDNase (100/500 reactions) DTT (Inactivation aid)
Quality control	Tested for absence of RNase
Storage	Store at ≤ -20°C

Protocol

The PRC decontamination protocol is designed for removing contaminating DNA from 20 µl reactions but can be scaled to other convenient volumes.



Contaminating bacterial DNA is present in most 2x qPCR master mixes

Contaminating bacterial DNA in commercial qPCR master mixes might cause false positive results when using qPCR for detection or quantification of minor amounts of bacterial DNA. We used an *E. coli* 23S primer/probe set to detect the presence of bacterial DNA in probe-based qPCR 2x master mixes from several suppliers. As shown in figure 1, contaminating bacterial DNA was found in all master mixes tested.

Decontamination of master mixes without reduction of sensitivity

Contaminating DNA is mainly a problem in high-sensitivity applications, where low-abundant DNA are targeted. Because of this, any loss of sensitivity in the qPCR assay caused by the decontamination protocol is unacceptable.




Fig 1. PCR master mixes are contaminated

PCR master mixes are contaminated. The presence of *E. coli* 23S DNA in 2x probe master mixes from various suppliers was quantified by using water as a no template control (NTC) and following the manufacturer's instructions. The figure shows plots acquired from several separate experiments. Traces of *E. coli* 23S DNA were found in all master mixes tested, with Cq values generally ranging from 30 – 35.




Fig 2. Decontamination with no impact on PCR sensitivity

Decontamination with no impact on PCR sensitivity. Untreated and decontaminated qPCR 2x master mix was used for analysis of an *E. coli* gDNA 10-fold serial dilution with 5 steps. NTC samples were included, and all plots of the serial dilution show an average of three replicates. Inset: Standard curves calculated from Cq values obtained from analysis of serial dilution.

In figure 2, the performance of decontaminated master mix was compared to untreated master mix.

A 5-step 10-fold serial dilution of *E. coli* gDNA was analysed using qPCR, which returned positive NTC. The PCR Decontamination Kit reduced the amount of contaminating DNA to levels below the detection limit of a 45-cycle qPCR experiment in the NTC samples. Treatment with the PCR Decontamination Kit did not significantly alter the Cq-levels obtained when analysing any of the dilutions of *E. coli* gDNA.

The Cq values presented was used to create standard curves (inset graph). No significant changes in either efficiency or R² between untreated and decontaminated master mix were observed.

Heat&Run® gDNA Removal Kit

- ✓ Single-tube protocol
- ✓ Suitable for qPCR settings
- ✓ Increase cDNA yield



“For removal of gDNA contamination from RNA preparations”

Genomic DNA can be a challenge when purifying RNA. Many RNA purification methods leave traces of gDNA, which can cause problems in later analysis.

Description

The Heat&Run gDNA Removal Kit is designed for removal of contaminating gDNA from high-quality RNA prior to reverse transcription in qPCR.

The kit is based on the recombinant heat labile dsDNase (HL-dsDNase), which is irreversibly inactivated at low temperatures. After DNase treatment the HL-dsDNase is easily inactivated by a moderate rise in temperature (58°C). Genomic DNA removal and reverse transcription can be performed in the same tube, thereby minimizing pipetting steps and reducing hands-on time.

Specifications

Source	HL-dsDNase is an engineered version of dsDNase
Contents	HL-dsDNase (50/250 reactions) 10X Reaction buffer
Quality control	Tested for absence of RNase
Storage	Store at ≤ -20°C


Protocol

Heat&Run gDNA Removal Kit can be used for RNA preps in samples with high-quality RNA, see figure below.

Purification of nucleic acids is required prior nuclease treatment.

Mix your RNA with Heat&Run, incubate for 10 min at 37°C, and then inactivate the dsDNase for 5 min at 58°C.

Run your RT!



No license required

At ArcticZymes, we pride ourselves on always offering seamless accessibility to our high-quality products. Produced under ISO 13485, our enzymes are sold under a “no license required” policy to ensure that our

customers are not restricted by legal burdens, now or with their future use. In addition, we offer our kits in a flexible format and are readily available to discuss your custom needs.

Ordering information

	Article No.	Pack size
PCR Decontamination Kit	80400-100	100 rxn
	80400-500	500 rxn
Heat&Run® gDNA Removal Kit	80200-50	50 rxn
	80200-250	250 rxn

Your OEM Partner to deliver novel solutions for genomics and proteomics

Quality

ArcticZymes is dedicated to the quality of our products and certified according to ISO 13485:2016.

Additional information

We are pleased to provide support using our protocols or customized protocols. For more information, please check our website www.arcticzymes.com or contact us.

Cutting-edge enzymes from Norway

ArcticZymes Technologies has a long history dating back to the late 1980s. Based in Tromsø, Northern Norway, we use access to the marine Arctic to identify new cold-adapted enzymes for use in molecular research, *in vitro* diagnostics and therapeutics. We focus on strong and reliable relationships with our business partners and commercial innovators around the world. Therefore, we are constantly striving to work at the highest level and not only meet but exceed the expectations of our partners.

In the service of science

The knowledge of the important role our enzymes play in research, diagnostics and therapeutics drives us every day. Our team of highly motivated and experienced scientists is constantly developing further innovations in order to expand our portfolio of novel and high-quality solutions.



A partner you can trust



Security of supply

With us you are always on the safe side when it comes to the timely delivery of high-quality enzymes. We strive for a reliable and uninterrupted supply of whatever enzyme technology you need.



Partnership approach

Our focus is on cooperative B2B partnerships which means that we put our customers' needs at the center of what we do. We strive to provide innovative solutions in order to help them to succeed in whatever they do.



Unique enzyme features

Enzymes play a decisive role in molecular research, *in vitro* diagnostics and therapeutics.

This makes it all the more important that they have a consistently high quality. Our novel enzymes are reproducible and have unique properties that make them particularly robust.



Unique access

Direct access to unique and diverse resources for bioprospecting allows us to continuously develop novel enzyme technologies with unique features and make them available to our partners.

Disclaimer: These products are intended for further manufacturing use or research use only. Certain applications of ArcticZymes Technologies ASA products may require licenses from others. It is the expressed duty of any receiver of ArcticZymes Technologies ASA products to acquire such licenses, if necessary. To the extent allowed by law, ArcticZymes Technologies ASA will not be liable for damages, whether direct, indirect, incidental, or consequential in connection with or arising from this document, including the use of it. ArcticZymes Technologies ASA products may be covered by pending or issued patents, designs or design applications and/or trademarks or trademark applications or any other registered or unregistered Intellectual Property Right. Version 1.0 • June 2021

ArcticZymes Technologies ASA

Sykehusveien 23
N-9294 Tromsø, Norway

T (47) 7764 8900
E contact@arcticzymes.com
I www.arcticzymes.com

ArcticZymes Technologies ASA

987 Old Eagle School Road, Suite 709
Wayne, PA 19087, USA

T (484) 534 3567
F (484) 368 3558
E contact-us@arcticzymes.com

ArcticZymes AS

Handelswei 1A 8501 XJ
Joure, the Netherlands

T (47) 7764 8900
E contact@arcticzymes.com
I www.arcticzymes.com

