

IsoPol™ SD⁺

IsoPol™ SD⁺ is the second in the ArcticZymes family of polymerases established to serve the genomic and proteomic market segment.

IsoPol™ SD⁺ is active at ambient temperatures. It has stronger strand displacement and higher salt tolerance, compared to similar polymerases and can be heat inactivated at temperatures above 50°C.

Features

IsoPol™ SD⁺ exhibits polymerase activity and lacks 3'-5' exonuclease activity. The 5'-3' exonuclease domain is truncated, thus the enzyme lacks 5'-3' exonuclease activity.

Optimal reaction conditions

IsoPol™ SD⁺ exhibits a specific activity of 7 500 U/mg. Optimal performance is achieved in a pH 8.5 buffer supplemented with 4-6 mM MgCl₂ and 100-350 mM salt (NaCl or KCl). Temperature optimum is at 37°C.

Inactivation

Inactivation is achieved by incubating the enzyme at 50°C for 20 minutes. It is also possible to inactivate the polymerase by addition of EDTA.

Unit definition

One unit is defined as the amount of enzyme that will incorporate 10 nmol of dNTP into acid insoluble material in 30 minutes at 37°C.

Strand-displacement activity

IsoPol™ SD⁺ exhibits improved strand-displacement activity compared to other similar DNA Polymerases at 25°C and 37°C.

Storage and stability

The enzyme is stable at -20°C for 2 years in the supplied storage buffer. Storage of the enzyme at 200 mM NaCl is recommended for high concentration samples.

Quality control

ssDNA endonuclease activity

25 U IsoPol™ SD⁺ was incubated with M13 ssDNA (0.5 µg) for 4 hours at 37°C. Agarose gel electrophoresis did not reveal any visible signs of ssDNA degradation.

dsDNA endonuclease activity

25 U IsoPol™ SD⁺ was incubated with a supercoiled plasmid (1 µg) for 4 hours at 37°C. Agarose gel electrophoresis did not reveal any transformation of closed circular DNA to nicked DNA.

ArcticZymes is dedicated to the quality of our products.

IsoPol™ SD⁺ is manufactured at our ISO 13485 certified facility in Norway.

Additional data and information

We are pleased to provide data and information relating to IsoPol™ SD⁺. Available data includes; buffer storage conditions, optimal salt, Mg²⁺, pH, processivity, activity, and strand displacement data.

Ordering information

Product name	Catalogue #	Concentration	Size	Units
IsoPol™ SD ⁺	71501-201	5 U/µl	40 µl	200 U
IsoPol™ SD ⁺	71501-100	According to agreement	According to agreement	According to agreement

Your OEM partner to deliver novel solutions for genomics and proteomics

ArcticZymes AS
Sykehusveien 23
N-9019 Tromsø, Norway

Phone: (47) 7764 8900
E-mail: contact@arcticzymes.com
Website: www.arcticzymes.com

ArcticZymes
600 West Germantown Pike, Suite 110
Plymouth Meeting, PA 19462, USA

Phone: (484) 534 3567
Fax: (484) 368 3558
E-mail: contact-us@arcticzymes.com

ArcticZymes
Handelswei 1A 8501 XJ,
Joure, the Netherlands

Phone: (47) 7764 8900
E-mail: contact@arcticzymes.com
Website: www.arcticzymes.com

Disclaimer: These products are intended for further manufacturing use or research use only. Certain applications of ArcticZymes AS products may require licenses from others. It is the expressed duty of any receiver of ArcticZymes AS products to acquire such licenses, if necessary. To the extent allowed by law, ArcticZymes AS will not be liable for damages, whether direct, indirect, incidental, or consequential in connection with or arising from this document, including the use of it. ArcticZymes AS products may be covered by pending or issued patents, designs or design applications and/or trademarks or trademark applications or any other registered or unregistered Intellectual Property Right.

Version 1.0 • December 2017